



ENDODONTIC EQUIPMENT

ENDO-MATE TC2

CORDLESS ENDODONTIC HANDPIECE
WITH TORQUE CONTROL & AUTO REVERSE

- ✦ Simple 5-key operation
- ✦ Light weight
- ✦ Cordless unit
- ✦ Alarm function alerts
- ✦ Variable torque
- ✦ Variable speed
- ✦ 5 programs
- ✦ Max. Torque : 4 Ncm (when using 20:1 head)
- ✦ Speed range: 100-2,250 min⁻¹ (With 20:1, 16:1, 10:1, 4:1 geared heads)
- ✦ Charge/continuous use: 90min/1.2h (depending on workload)

