

- ✦ Independent up & down and backrest movement of the chair
- ✦ Chair interlock system, synchronized chair
- ✦ Big upholstery with chair interlock system, more safe and comfortable
- ✦ 9 customised programs for 3x3 users
- ✦ Luxury 6 LED sensor based light 10,000 - 25,000 Lux
- ✦ Safety switch for water - air - electricity with single on/off switch control
- ✦ Built in type chair (without floor box)

**DELIVERY UNIT**

- ✦ With pneumatic locking arm to prevent any unnecessary displacement of the unit
- ✦ Feather touch control panel for easy touch play operation
- ✦ 3 air rotor points and extra empty/ vacant point
- ✦ One 3 Way syringe
- ✦ Detachable and easy clean handpiece holder block
- ✦ Control panel with 16 keys for chair control and other functions
- ✦ LED X-Ray viewer
- ✦ LCD display
- ✦ Inbuilt woodpecker scaler

**WATER UNIT**

- ✦ Rotatable water unit & spittoon
- ✦ Glass cuspidor
- ✦ Detachable cuspidor, easy to clean
- ✦ Movable water unit for more space to the assistant
- ✦ Built out and detachable suction filter for easy to clean
- ✦ Tissue holder
- ✦ Water heater for lukewarm water



**GALAXY+**



Rotatable Utility Tray

EN ISO 13485:2012



**MULTIFUNCTION FOOT CONTROL**

- ✦ Unique metal handle
- ✦ Chair position control
- ✦ Bowl rising control
- ✦ Cup filling control
- ✦ 2 independent switches for air and water
- ✦ Stainless steel tubing for foot control operation

**ASSISTANT ARM**

- ✦ Low vacuum and high vacuum cannulas
- ✦ 3 Way syringe with option for Luke warm water
- ✦ Chair control for assistant

**FULL FUNCTION CONTROL PANEL WITH LCD SCREEN & 16 BUTTONS**

- Setting
- Cuspidor flush
- Cup filling
- Water heating
- Operation light
- Zero position
- Spitting position
- Doctor's memory selected
- Chair position 3 memory for 1st user
- Chair position 3 memory for 2nd user
- Chair position 3 memory for 3rd user
- X-Ray film viewer

**SPITTING POSITION SETTING**

On pressing the spitting position button, automatically the backrest moves upward, bowl rinsing start, light turns off, if pressed again, chair would automatically return to working status, backrest backward, bowl rinsing closes, lights turns on

