



# Surgic Pro<sup>+</sup>


## Optic

- ✦ Surgic Pro+ with 230W high power
- ✦ 5~80 Ncm powerful torque
- ✦ Wide speed range of 200 ~ 40,000 min<sup>-1</sup>
- ✦ Large backlit LCD Display to see & check the displayed data
- ✦ User can save up to 8 preferred settings (speed, torque, coolant flow and rotational direction)
- ✦ Control Unit with advanced hand piece calibration for torque & speed accuracy
- ✦ **Data Accessibility** - Treatment data can be accessed & downloaded using a USB memory stick. (USB memory stick is **optional**)
- ✦ Micro motor with 2 meter Cord
- ✦ **Data storage facility** - Can record the date, speed, rotational direction, torque, gear ratio, coolant solution flow volume up to 100min, with no special software needed.
- ✦ With X-DSG20L dismantable Hand piece (Optic)
- ✦ Control panel with hand piece light Intensity variation control (Off / Medium / High)
- ✦ Light power : Over 32,000 LUX (LED Micromotor)
- ✦ Dimensions : W268 x D220 x H100 (mm)
- ✦ Built in memory : 100 min

## IMPLANT MICROMOTOR



**Control Unit with  
Advanced Handpiece  
Calibration (AHC) -  
High Torque Accuracy**



### Foot Control :

- ✦ Speed control
- ✦ Coolant control
- ✦ Program control
- ✦ Forward/Reverse direction
- ✦ Speed control : Progressive



**Simple and Intuitive  
Operation Control Panel**

