



## NSK HANDPIECE

# DynaLED

Self Illuminating LED

- ✦ Connect to any standard air tubing coupling with any dental chair
- ✦ **No need of circuit board, tubings & bulbs**
- ✦ High intensity & stable self LED illumination through mini power generator
- ✦ Cellular glass optic - pure daylight quality illumination
- ✦ Stainless steel : durable, stable & light weight
- ✦ Available with Miniature Head M500LG (380,000 - 450,000 RPM) Head Size:  $\varnothing$ 10.6 x H12.4mm, 25,000 LUX & Standard Head M600LG (330,000 - 430,000 RPM) Head Size:  $\varnothing$ 12.1 x H13.3mm, 25,000 LUX
- ✦ **Clean head system**
- ✦ **Push button chuck**
- ✦ **Quattro Spray**

### Visibility & Access :

Improved visibility and direct access to the treatment field allows faster and more accurate preparations.

### LED Durability :

The LED light source is significantly more durable than quartz. Replicated daylight illumination quality.



Optic - NSK LED



Non Optic