

FLARE DENTAL CHAIR

- ✦ Independent up & down movement and backrest movement of the chair
- ✦ Luxury LED operation lamp
- ✦ Right arm movable for easy seating of the patient
- ✦ 2 programs in chair ...1 preset working position other for Zero position
- ✦ Feather touch control panel present on delivery unit
- ✦ Doctor's stool with Adjustable height

MOVEABLE FOOT CONTROL

- ✦ 1 Pedal for chip air (Blow)
- ✦ Pedal for air
- ✦ 1 Pedal for air & water
- ✦ Round joystick on chair base for chair movements

ASSISTANT ARM

- ✦ Movable assistant arm
- ✦ 2 suction point for low and medium suction
- ✦ 1 - 3 way syringe



CE ISO 13485 (2003) / AC:2007 / AC:2009



TRADITIONAL DELIVERY UNIT

- ✦ Pneumatic locking Arm
- ✦ Feather touch control panel for all controls of chair
- ✦ 1 - 3 way syringe
- ✦ 3 points for air rotor
- ✦ 1 extra holder
- ✦ A pressure gauze present for the air rotor
- ✦ 1 LED X-Ray viewer

WATER UNIT

- ✦ Fixed water unit with movable glass spittoon
- ✦ Cup filler with option for warm & normal water