



DIGITAL PIEZON ULTRASONIC SCALER

D7 LED SCALER

- ✦ Detachable handpiece with light, better visibility, autoclavable
- ✦ Works on 3 modes : Scaling, Perio & Endo
- ✦ Clinical solutions are applicable for automatic water supply mode, improving the performance of perio and endo function dramatically
- ✦ Less heat from tip during operation. Improve comfort degree of patients
- ✦ Output power 3W~20W
- ✦ Frequency 28kHz \pm 3kHz
- ✦ Dimensions 237mm x 206mm x 95mm



Fiber optic (LED) Handpiece