

VS Systems

Hygienic & Ergonomic Access to preparation

- ✦ Continuous high suction power
- ✦ Robust VS suction technology with automatic separator to separate liquid from solids
- ✦ Guarantees high reliability
- ✦ High suction performance within seconds
- ✦ Flexible installation by custom-made program & functional accessories



VS 250 S



VS 300 S



VS 600

Data & Facts			
Model	VS 250 S	VS 300 S	VS 600
Voltage (V)	230 (1~)	230	230/400
Frequency (Hz)	50 / 60	50/60	50/ 50/60
Current consumption (A)	90 at 250 l/min	2.9/3.7	5.0/ 1.8 - 4.0
No. of treatment chairs	1	-	
Power (kW)	-	0.58/0.8	1.1/ 1.1/1.4
Speed (rpm)	-	2,850	2,850
Vacuum (mbar)	-	180	
Cyclic duration factor (S1)	-	100%	100 %
Max fluid flow rate (l/min)	-	4	-
Noise level without housing dB(A)	64 / 66	Approx. 63-65	Approx. 63/ Approx. 63
Noise level with housing dB(A)	52 / 54	Approx. 54-56	Approx. 55/ Approx. 55
Dimensions without housing (H x W x D cm)	32 x 28 x 28	37 x 31 x 31	50 x 42 x 41
Dimensions with housing (H x W x D cm)	47 x 35 x 46	51 x 36 x 36	65 x 49 x 53
Weight without housing (kg)	-	12.5	25
Weight with housing (kg)	-	20	45

 Made in Germany