



## DIGITAL PIEZON ULTRASONIC SCALER

### DTE - D5

#### Optic / Non-Optic

- ✦ Detachable & autoclavable handpiece
- ✦ Works on 3 modes : General, Perio & Endo
- ✦ Automatic frequency tracking
- ✦ With 5 Titanium Scaling tips (3 S, 1 P, 1 Stain Removal) & 1 Endo attachment
- ✦ Power input : 100V to 240V~ 50Hz/60Hz 1.2A (Max)
- ✦ Output Power : 3W to 20W
- ✦ Frequency : 28 kHz  $\pm$  3 kHz
- ✦ Dimension : 188 mm X 159 mm X 81 mm
- ✦ Main unit weight : 0.65 kg



Fiber optic (LED) Handpiece