

Variable Frequency Suction Machine

- ✦ Wet line suction system
- ✦ High suction power, stable & concentrated
- ✦ Internal noise reduction system
- ✦ Use of frequency conversion technology, can adjust the flow according to the chairs that are actually being used



Data & Facts			
Unit Type	VMAX 450	VMAX 1250	VMAX 1450
Power (HP)	1	2	3
Max Pumping (L/min)	500	1000	1500
Vacuum Degree (KPa)	-13	-21 ~ -24	-22 ~ -25
Suitable Chairs	Upto 2 Chair	Upto 4 Chair	Upto 6 Chair
Noise (dB)	53	57~62	58~62