



DIGITAL PIEZON ULTRASONIC SCALER

DTE - D3

Optic / Non-Optic

- ✦ Detachable & autoclavable handpiece
- ✦ Works on 2 modes : General & Perio
- ✦ Automatic frequency tracking
- ✦ With 5 Tips (3 S, 1 P, 1 Stain Removal)
- ✦ Power input : 100V to 240V~ 50Hz/60Hz 1.2A (Max)
- ✦ Output power : 3W to 20W
- ✦ Frequency : 28kHz \pm 3kHz
- ✦ Dimension : 196 mm X 134 mm X 80 mm
- ✦ Main unit weight : 0.62 kg



Fiber optic (LED) Handpiece