



DIGITAL PIEZON ULTRASONIC SCALER

INBUILT SCALER V2

- ✦ Detachable handpiece
- ✦ 2 modes
- ✦ With 5 Tips (3 S, 2 P)
- ✦ With Torque wrench
- ✦ Power input : 100V to 240V~ 50Hz/60Hz
1.2A(Max)
- ✦ Output Power : 3W to 20W
- ✦ Frequency : 28kHz \pm 3kHz
- ✦ Dimension : 74 mm X 56 mm X 38 mm
- ✦ Main unit weight : 0.2 kg
- ✦ Inbuilt Scaler V1 LED & V3 LED also available

