



LIGHT CURE UNIT

LIGHT CURE UNIT

- ✦ Strong optical output power
- ✦ Three treatment modes : Continuous, Ramp & Pulse
- ✦ Digital time setting : 0~20s for LED curing light
- ✦ Full color LCD display
- ✦ LED cold light output can treat all dental tissue
- ✦ Alarm for low voltage

Technical Specifications

- ✦ Working voltage : 100-220 V AC, 50-60 Hz
- ✦ Working volt : 4.2-4.5 volts
- ✦ Electric current : 0.7 amperes
- ✦ Luminous intensity : > 1200 mw/cm² 5W
- ✦ Wavelength : 470 nm
- ✦ Weight : 160g
- ✦ With 1 black probe



EC ISO 13485:2012 / AC:2012