



ADVANCE PIEZOSURGERY UNIT

HIGH-PERFORMANCE HANDPIECE

Optimized handpiece with stronger sealing technology. High temperature cooking disinfection causes no loss of performance. High brightness LED brings clearer operating vision

Clinical applications

Minimally invasive tooth extraction, Intralift, Sinus Lift, Alveoplasty, apical cyst extraction, Crown Extension, Crest splitting etc.

Clinical Benefits

- ✦ Micron Cutting
- ✦ Selective Cutting With No Injury to Soft Tissue
- ✦ Cold Cutting Mode

For different treatment needs, Four different working modes are set :

- ✦ Cutting Mode
- ✦ Periodontal Mode
- ✦ Endodontic Mode
- ✦ Cleaning Mode

Technical Specification :

- ✦ Dimension : 276mm x 267mm x 110mm
- ✦ Intermittent running device : Intermittent running with 10s pause after 60s operation
- ✦ Water flow : 25~110ml/m
- ✦ Mains Voltage : 100-240V~50Hz/60Hz 150VA
- ✦ Working Frequency : 24-36KHZ
- ✦ Weight : 2.8KG

SURGERY X

The Third Generation Woodpecker Ultrasurgery



70
Watts

5 Yrs
warranty
on PCB

EQUIPPED WITH 2 FIBER OPTIC HANDPIECES & 14 TIPS (BASIC KIT + SINUS LIFT KIT)

NEW ADVANCE FEATURES

Tips Kit

Specially used for tooth extraction, crest splitting & retro preparation. Adopt gold plating which is wear-resistant & practical with long time of endurance.



7-Inch touch screen

Abandoning the previous button operation settings, innovatively adopt 7-inch screen. While clearly displaying information, operate directly by tapping the screen.



Multi-function Foot Pedal

Foot Pedal makes it easy to control water flow, working mode and power, simplifying operation procedures and realizing efficient operation.



14 Tips

US5 US4 US3 US2 US1R US1L US1 UL5 UL4 UL3 UL2 UL1 UC1 UI2



Dental Chair
Compressor & Suction
Autoclave & Piezosurgery Unit
Endomotor & Apex Locator
Implant Micromotor & Piezosurgery Unit
Clinical Equipment