

- ✦ Simple, easy & quick to use
- ✦ Convenient cordless device using battery
- ✦ Cordless & electric motorized filling system
- ✦ Refill Cartridge
- ✦ 360° rotation at both cartridge needle & pack tip
- ✦ High capacity lithium-ion battery

Clinical Benefit

- ✦ Less time consuming
- ✦ Less microbial leakage
- ✦ Adapts better to grooves & depressions of coronal Walls & Lateral Canals
- ✦ No dead space preventing secondary infection

Pack :

- ✦ Instant heating & instant cooling of Pack tip
- ✦ 360° rotation without taking out the pack tip
- ✦ Outstanding heating performance (Max Temp. : 300°C)
- ✦ Ergonomically designed body for comfortable grip
- ✦ Highly visible LED display
- ✦ "Hold" switch for safety
- ✦ Auto off : 5min
- ✦ Dimension : 153mm x 27mm x 25mm (WxLxH)
- ✦ Weight : 90 g (Including battery)
- ✦ Operating Voltage : 3.6V
- ✦ Range of setting temperature : 140, 200, 300 °C

Fill :

- ✦ Easy-Install refill needle
- ✦ High visible LED display
- ✦ 360° rotating cartridge needle
- ✦ Fully motorized obturation for smooth & easy filling
- ✦ Ergonomically designed body for comfortable grip
- ✦ Auto off : 5min
- ✦ Auto reverse function when the cartridge is empty
- ✦ "Hold" switch for safety
- ✦ Dimension : 223mm x 27mm x 26mm (WxLxH)
- ✦ Weight : 177 g (Including battery)
- ✦ Operating Voltage : 3.6V
- ✦ Range of setting temperature : 140,160,180, 200 °C

Genuine System, Genius Result

GENESYS

**Genuine Root Canal
Obturation system
for Thermal 3D
Obturation of Root Canal**



**Pack handpiece for
Down Packing**



**Filler handpiece for
Back Fill**

