



LIGHT CURE UNIT

O-LIGHT LIGHT CURE UNIT

- ✦ **1 Sec. for curing** - Cure 2mm resin in 1 second, with maximum light intensity above 2300mw/cm², 5 W high-power blue light LED guarantees light intensity output
- ✦ **High power LED lights** - Imported wick from America, 30% lower current under same light intensity. Lower heating, longer service life.
- ✦ **Ergonomic, feel great in hand** - Smooth body line, Good hand-holding feeling, clear shape lines, classical back & white color matching
- ✦ **Strong penetration** - Strong penetration, without distance limitation. Long distance irradiation does not affect Curing effect
- ✦ **Flexibility** - 360° rotation of front connector, easy to reach areas where used to be hard to reach. Enable flexible operation in posterior teeth area. Small mouth-open during operation brings patients comfortable treatment. 90° vertically irradiate operation area bringing higher curing efficiency
- ✦ **Two modes for option** - High Mode is of 2300mw/cm²- 2500mw/cm² light intensity is more suitable for orthodontics, porcelain veneer & post-core restoration. Standard Mode is of 1000mw/cm²- 1200mw/cm² is suitable for general curing

