

## LIGHT CURE UNIT

### DIGITAL LED

### LIGHT CURE UNIT

- ✦ Material of shell handle : Aluminium
  - ✦ Size : Diameter 23 X 160mm
  - ✦ Weight : 160g
  - ✦ Grade : Ordinary B Type Intermittent working
  - ✦ Working Voltage : 100-220V AC, 50-60HZ
  - ✦ Grade : II IP40
  - ✦ Mode : Three kinds of working mode (Alterability)
  - ✦ Working Volt : 4.2~4.5 volts
  - ✦ Electric current : 0.7 amperes
  - ✦ Battery Type : Lithium Ion
  - ✦ Battery capacity : 1400Mah
  - ✦ Wave Length : 420-480nm
- Luminous Intensity : 1200mW/cm<sup>2</sup> 5W

(According to IEC60601-1, curing lights belongs to Class II, I P40)



**TULIP**