

INTRA ORAL CAMERA WITH MONITOR

- ✦ Camera features a 1/4" High Definition image sensor
- ✦ 17" Monitor features a High Definition resolution of 1280 X 1024
- ✦ 8 GB Memory Card (included) allows for easy image transfer & file management
- ✦ Camera buttons are easily accessible to control capture & playback functions
- ✦ Camera is illuminated with six bright LEDs, giving the best contrast, definition & resolution quality
- ✦ Playback format includes : MP3 Audio, MP4 Video & JPG images
- ✦ Adjustments on monitor include : Contrast, brightness & color
- ✦ Monitor comes equipped with a medical grade tempered glass front panel to minimize scratches
- ✦ Easy integration of LCD monitor with camera on the dental delivery unit due to no additional cable needed
- ✦ Input supports USB, VGA, HDMI & AV

