

INTRA ORAL CAMERA

ADVANCE CAM IOC

VGA + USB + AV MODEL

- ✦ Intra Oral camera with 6 LED
- ✦ CCD Sensor
- ✦ Memory - Camera 50 pic & USB 10,000 pic
- ✦ VGA & USB 2.0 signal Output
- ✦ Resolution : 2 mega pixels
- ✦ Compatible to USB, TFT Monitor without CPU
- ✦ Lens 0.02 inch pin lens
- ✦ Focus Range 8-12 mm
- ✦ Backup - USB device backup
- ✦ Display modes - Full screen & 4 split modes
- ✦ Intelligent function - Prevent dithering, Anti-alias edge
- ✦ Now available with Window 7 (Compatible)

