

## INTRA ORAL CAMERA

### WIFI INTRA ORAL CAMERA

- ✦ End-user : Support 4 users at the same time
- ✦ Available Distance :  $\geq 30\text{m}$ (open space)
- ✦ Battery capacity : 2100mA/h
- ✦ Working time : 140mins
- ✦ Pixel : 1 million
- ✦ Light Resource : 6pcs white LED light (5600K)
- ✦ Dimming : Five brightness adjustment
- ✦ Focus Range : 10mm - 50mm
- ✦ Photography Angle :  $105^\circ$
- ✦ Connection : WiFi
- ✦ Materials : Imported medical resin
- ✦ Main Function : Real-time watching, oral examination
- ✦ Two-way take photo and video
- ✦ Compatible with iphone & Android devices



HD : 720P