

## CLINICAL EQUIPMENT

### X-RAY FILM VIEWER

- ✦ LED light source, working life over 100000 hours
- ✦ Can mount on the wall or place on the desk with 3 positions
- ✦ Fashionable design
- ✦ Intensity of light : >3000 LUX
- ✦ Dimensions : X-Ray viewer - 35,7 x 24,3 x 3,5 cm (W x H x D) / Screen - 30 x 20 cm (W x H)
- ✦ Illuminance : 7.000 Lux
- ✦ Light : Cold, polarized, no oscillation
- ✦ Power : 10 Watts
- ✦ Voltage : 220 V - Automatic voltage selection
- ✦ Net Weight : 1,1 kg



Cristófoli

Luna