

DENTAL CABINET

CLINICAL EQUIPMENT

MOBILE DENTAL CABINET

- ✦ High quality stainless steel with anti-fingerprint
- ✦ Different size drawers according to usage
- ✦ Concealed resistance track
- ✦ Artificial marble tabletop with wear & corrosion resistant, anti-penetration, no stain & easy to clean
- ✦ **Wheels with soft movements**
- ✦ **Stylish & durable handle**
- ✦ **Contains 5 different size - 5 nos. drawers**



Model No. Size in mm	Height	Width / Depth	Length
Mobile Dental Cabinet	820	500	500