

ULTRASONIC CLEANER

- + Tank Material : SUS304 Stainless Steel
- + Tank Capacity : 10 Ltr
- + Time setting :1-60 min adjustable
- + Temperature setting : 200° C - 800° C adjustable
- + High-performance transducer
- + CE/FCC/RoHS tested and approval
- + More quiet to operate in the hospital
- + Special groove design for better waterproof

DIGITAL ULTRASONIC CLEANER

- + 3 color LED display
- + 1 to 30 - minute full range timer
- + 5 temperature settings
- + Solution usage timer
- + Degas feature
- + Independent power switch and detachable power cord, convenient to use
- + Tank capacity 3 Ltr
- + Frequency 35,000 Hz
- + Tank Size 25 x 17.4 x 7.8 cm and Weight 5.1 Kg
- + Tank material stainless steel SUS304

MIXER- AMALGAMATOR

- + Digital device with ergonomic design, safe and highly Economical
- + Plastic armature and wide window
- + Digital display to set the mixing frequency and time selection (1-99s)
- + Uses pre-dose capsules of alloy and mercury
- + Time and oscillations can be varied with an amplitude 25nm & frequency of 4350 oscillations/ second
- + Easy to use, stable, extremely quiet, easy to clean and maintain

GT SONIC®



CE GB/T/9001-2008 ISO 9001:2008

CODYSON



Standard Accessories :



Basket

Optional Accessories :

- + SS Tray + Rack
- + Beaker + Beaker Holder

