



# Implant-X

Precision Strength Safety

Your favorite assistant of high-precision dental implanting treatment



Strategy Partner

**UNICORN DenMart**   
Behind The Best Indian Dentists

## Features

### ✦ Perfect Balance

Balance between weight of contra angle & motor guarantees & reduces fatigue in complex surgeries



### ✦ Outstanding Quality Makes It Your Solid Choice

Implant X adopts noiseless motor imported from Swiss



### ✦ Advance FOC (Field of Current)

Blends low speed performance with high speed performance

### ✦ Strong & Stable Contra-Angle

High precision LED handpiece with less vibrations enables stable torque output under high & low speed



### ✦ LED Contra Angle

Ensures visibility, high accuracy & efficiency

### ✦ GRW High Precision Miniature & Handpiece Bearings

Guarantees performance & durability of contra angle

### ✦ Color Touch Screen With Glass Surface

Easy to clean with disinfectants ensuring safety



### ✦ Intuitive User Interface

Large color touch screen displays every procedures in Implant surgery with pre set parameters eliminating setting process



### ✦ Mute Peristaltic Pump

Easy to connect & compatible with disposable water tubings

### ✦ Flexible Multi-Function Foot Pedal

For aseptic operations



## Specification

Rotation Speed Range	300-40,000 rpm
Torque Range/ Gear 20:1	5Ncm-80Ncm
Max Pump Output	110ml/min
Main Unit Dimensions	276mm*267mm*110mm
Motor Dimensions	φ 23.4mm*L91mm
Motor Weight	140g
Power Supply	AC220V 50/60Hz

Strategy Partner

**UNICORN**  
**DenMart**   
Behind The Best Indian Dentists