

PERFECT COMBINATION OF RELIABILITY & POWERFUL INNOVATION

PORT-X II

- ✦ Portable, compact & wireless DC X-Ray
- ✦ 60kV/2mA Toshiba Tube
- ✦ Focal Spot 0.8 mm
- ✦ Exposure Time (0.01-2.0 sec)
- ✦ Graphic LCD display
- ✦ IOPA film & Radiovisiography sensor compatible
- ✦ Round type beam limiting device
- ✦ Total filtration 1.8 mmAl
- ✦ 70-90 X-Ray shots on film or 100 to 120 X-Ray shots on Sensor once battery is recharged
- ✦ Weight approx. 2.8 Kg
- ✦ Electrical Voltage 220 V
- ✦ Internally LEAD coated to prevent excessive exposure & dispersion of radiation
- ✦ Rechargeable Lithium-Polymer battery
- ✦ Patient & tooth selection mode
- ✦ Very useful for multi-chair clinics

Compact Outside - Smart Inside



Digital Sensor 0.15 sec
Chemical Film 0.3 sec



Canine

Digital Sensor 0.15 sec
Chemical Film 0.3 sec



Incisor

Digital Sensor 0.25 sec
Chemical Film 0.35 sec



Molar

India's Highest selling
Portable X-Ray with
7000+ Installations