

Dental Sintering Furnace

KDF

Zircon

Made in Japan, the compact & user friendly Zircon Furnace from KDF is designed to sinter a variety of Zirconia materials.

With Quick Firing Program



- ✦ It features molybdenum disilicide heating elements which provide clean constant & reliable results.
- ✦ It can be programmed with 100 programs with 20 steps each.
- ✦ It has a quick firing program with a cycle of upto 1.5 hours with high temperature.
- ✦ It can sinter 75 units simultaneously with 3 stackable trays.
- ✦ Comes with SD cards which allows users to edit their firing programs

Technical Parameters

Model	KJ-M1700-DL
Firing Platform	Ø95mm (3.74 inch)
Working temperature	1600 °C
Display	3.5inch LCD Color Display
Heating Element	High Purity MoSi2
Working power	AC 220V. 50/60Hz. Max. 2.5kVA
Dimensions	12 (w) x 25 (h) x 18 (d) inches
Gross weight	40Kg